DXR RAMAN INSTRUMENTS

Safety Guide

The information in this publication is provided for reference only. All information contained in this publication is believed to be correct and complete. Thermo Fisher Scientific shall not be liable for errors contained herein nor for incidental or consequential damages in connection with the furnishing, performance or use of this material. All product specifications, as well as the information contained in this publication, are subject to change without notice.

This publication may contain or reference information and products protected by copyrights or patents and does not convey any license under our patent rights, nor the rights of others. We do not assume any liability arising out of any infringements of patents or other rights of third parties.

We make no warranty of any kind with regard to this material, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Customers are ultimately responsible for validation of their systems.

© 2008 Thermo Fisher Scientific Inc. All rights reserved. No part of this publication may be stored in a retrieval system, transmitted, or reproduced in any way, including but not limited to photocopy, photograph, magnetic or other record, without our prior written permission.

For technical assistance, please contact:

Technical Support
Thermo Fisher Scientific
5225 Verona Road
Madison, WI 53711-4495
U.S.A.

Telephone: 1 800 532 4752 (U.S.A.) or +1 608 273 5017 (worldwide)

Fax: +1 608 273 5045 (worldwide)

E-mail: us.techsupport.analyze@thermofisher.com (U.S.A.)

World Wide Web: <http://www.thermo.com/spectroscopy>

269-215200, Rev A

Contents

Introduction	1
Questions or concerns.....	1
Manual conventions.....	1
General safety.....	3
Safety labeling.....	3
Fire safety and burn hazards	3
Environmental conditions.....	4
Lifting your instrument.....	4
Corrosives and solvents	4
WEEE compliance.....	4
Electrical Safety.....	5
Electrical service.....	5
Fuses	6
Illuminator bulbs	6
Grounding.....	6
Power cords	7
Laser Safety	7
Laser safety precautions.....	8
Laser emissions.....	9
Manufacturer's laser information	10
Laser safety goggles.....	11
Disabling the excitation laser.....	11
Protective housing.....	11

Introduction

This manual contains the safety precautions that must be followed when using a Thermo Scientific DXR Raman instrument. These instruments are intended to be used to analyze the chemical composition of materials using dispersive Raman spectroscopy. Using these instruments in a manner not specified in the documentation that came with the instrument may create a hazardous condition. Each person who will be using the instrument should read this manual. Safety information is also included in the printed and on-line documentation that came with the instrument.

Questions or concerns

In case of emergency, follow the procedures established by your facility. If you have questions or concerns about safety or need assistance with operation, repairs or replacement parts, you can contact our sales or service representative in your area or use the information at the beginning of this document to contact us.

Conventions

This manual includes safety precautions and other important information presented in the following format:

- Note** Notes contain helpful supplementary information. ▲
- Notice** Follow instructions labeled “Notice” to avoid damaging the system hardware or losing data. ▲
- ⚠ Caution** Indicates a hazardous situation which, if not avoided, could result in minor or moderate injury. ▲
- ⚠ Warning** Indicates a hazardous situation which, if not avoided, could result in death or serious injury. ▲
- ⚠ Danger** Indicates a hazardous situation which, if not avoided, will result in death or serious injury. ▲

The following safety symbols are used on this product:

<i>Symbol</i>	<i>Description</i>	<i>Indication</i>
	Black graphical symbol inside a yellow triangle with a black triangular band	This is warning symbol. The graphic in this symbol is used to alert the user to potential hazards.
	Black graphical symbol inside a red circular band with a red diagonal bar	This is a prohibition symbol. The graphic in this symbol is used to alert the user to actions that shall not be taken or shall be stopped.
	White graphical symbol inside a blue circle	This is a mandatory action symbol. It is used to indicate that an action shall be taken to avoid a hazard.
	Black graphical symbol inside a yellow triangle with a black triangular band	This is the general warning sign. Failure to heed the safety precautions could result in personal injury.
	White graphical symbol inside a blue circle	This is the general data loss or property damage symbol and is not related to personal injury. Failure to heed these precautions can result in irreparable damage to property or permanent data loss.

General safety

The following sections contain general safety information and advice you should always follow when using your instrument.

Caution

When operating your instrument, use only accessories and equipment supplied by us. Use of equipment and accessories not supplied by us can result in permanent damage to the instrument and may create a safety hazard. ▲

Safety labeling

In many cases, safety information is displayed on the instrument itself. If you see the symbol shown here on an instrument or accessory, be sure to read the safety label and check the documentation that came with the instrument or accessory. The symbol shown here indicates that there is additional information in the documentation. Additional safety labels are shown throughout this manual and in the on-line instrument help system. If any safety labels become loose or unreadable, contact us to obtain replacement labels.

Fire safety and burn hazards

The following symbols appear inside the spectrometer to warn you about hot surfaces that could cause a burn injury in the vicinity of the symbol.

Caution

To avoid the risk of fire or explosion, do not test flammable or explosive samples. ▲

Caution

To avoid the risk of fire, never block any of the vents on the instrument or the power supply. ▲

Environmental conditions

Our instruments and accessories are designed for indoor use at altitudes up to 2,000 meters (6,500 feet). They operate reliably at temperatures between 16 and 27 degrees C (60 and 80 degrees F). Maintain humidity in the range of 20% to 80%.

Lifting your instrument

Your instrument weighs between 46-57 kg (100-125 Lb). To avoid risk of injury, one person should not attempt to lift the instrument. It should be lifted by at least two people, and safe lifting techniques should always be used.

Corrosives and solvents

Many standard spectroscopy methods are based on the use of solvents. Sample materials dissolved in solvents can be measured using your spectrometer, but special precautions must be taken.

▲ Warning

Halogenated hydrocarbons can oxidize when exposed to high temperatures in the instrument. This can produce materials such as hydrochloric acid, hydrofluoric acid and phosgene, which are highly toxic and corrosive. Materials like these can cause personal injury or death and can damage the instrument. If you regularly use solvents containing halogenated hydrocarbons, be sure your work area is properly ventilated. ▲

WEEE compliance

Many instruments are required to comply with the European Union's Waste Electrical & Electronic Equipment (WEEE) Directive 2002/96/EC. If compliance is required, the instrument is marked with the symbol shown here.

We have contracted with one or more recycling/disposal companies in each EU Member State, and this product should be disposed of or recycled through them. Further information is available at www.thermo.com/WEEERoHS.

Electrical Safety

The following sections contain electrical safety information and advice you should always follow when using your instrument. If you see the symbol shown here, there is a risk of electric shock in the vicinity of the symbol. Only qualified service representatives from our company should attempt to service a component that carries this symbol.

⚠ Warning

Avoid shock hazard. Do not allow liquid to run into the power supply. ▲

⚠ Warning

If a protective cover on the instrument, computer, or monitor appears damaged, turn off the system and secure it against any operation. Always examine any protective covers for transport stresses after shipping. ▲

⚠ Warning

Even after the instrument has been disconnected from all voltage sources, capacitors may remain charged for up to 30 seconds and can cause an electrical shock. ▲

⚠ Caution

Do not position your instrument in a way that would make the power switch difficult to access or operate. ▲

Electrical service

Check the current, voltage, and line frequency of the wall outlet that will be used for the instrument *before* you connect the instrument or turn on the power. (This should be done before or during the system installation and any time the system is plugged into a different wall outlet.) In addition, your electrical service must also be stable. The following table provides the specifications for line disturbances and noise, as well as current, voltage, and frequency information.

⚠ Danger

To avoid injury, only a qualified person using the appropriate measuring device should check the line voltage, current and line frequency. ▲

<i>Characteristic</i>	<i>Specification</i>
Input current	1.7 A at 230 VAC
Input voltage	90-264 VAC
Line frequency	47-63 Hz
Line disturbances	Sags, surges, or other line disturbances must not exceed 10% of input voltage (even for a half cycle)
Noise	less than 2 V (common mode) less than 20 V (normal mode)

If you are not sure that your electrical service meets these requirements, contact us for information about a power line audit.

Fuses

The fuses for the spectrometer and microscope are in the power supply and cannot be replaced. If you suspect that a fuse has blown, contact your service representative for assistance.

Illuminator bulbs

The standard microscope illuminator bulb is an LED that should not burn out during the lifetime of the microscope. If it does, or if there is a problem with the illuminator or illuminator bulb, contact your service representative for assistance. If you have an optional illuminator with a bulb that can be replaced, refer to the manufacturer's documentation for instructions for replacing the bulb.

Grounding

Each wall outlet you use must be equipped with a 3-wire line: live, neutral, and ground. The ground must be a non-current-carrying wire connected to earth ground at the main distribution box. To assure a good ground connection and avoid shock hazard, do not use an outlet that has ground connected to a conduit ground.

Power cords

You must be sure to use an appropriate power cord for the electrical service. The power cord supplied with the spectrometer is a 3-wire, grounded power cord, appropriate for use in the country listed as the shipping destination for the spectrometer. If the power cord you received is not appropriate for the electrical system in your location, contact us to order a new power cord.

If the power cord becomes damaged, replace it. Contact us if you need additional information about replacement cord or extension cord specifications.

Warning

To prevent electrical hazards, do not remove or defeat the ground prong on the power cord, and do not use any cords that are frayed or damaged. If you use an extension cord, it also must have a protective conductor. ▲

Laser Safety

Your instrument may be a class I or class IIIb laser product, depending on the configuration of the instrument. Exposure to any laser radiation can be dangerous and should be avoided. Class IIIb laser radiation is particularly dangerous and can cause severe burns and permanent eye damage. Anyone working with or around one of our Raman instruments (especially a class IIIb configuration) should follow all national or local laser safety requirements. All such personnel should also read and follow the laser safety information and warnings in this chapter and should read the “Setting up a laser safety area” topic in the on-line instrument help system to help avoid the risk of injury. If you see the symbol shown here, there is a risk of exposure to visible or invisible laser radiation in the vicinity.

- ⚠ Warning** Verify that your laboratory is compliant with the American National Standard for the safe use of lasers (ANSI Z136.1-1986) and any other national or local laser safety requirements. This includes, but is not limited to:
- Hazard evaluation and classification
 - Control measures
 - Laser safety and training programs
 - Criteria for exposure of the eye and skin. ▲

Laser safety precautions

The United States Department of Health and Human Services warns against improper laser use, as follows:

- ⚠ Warning** A class IIIb laser product should always be operated in a laser safety control area. (See the “Setting up a laser safety area” topic in the on-line help spectrometer system for more information.) All persons using or in the vicinity of a class IIIb laser product should be aware of the potential hazards and should wear appropriate laser safety eyewear. (See “Laser safety goggles” in this chapter for more information.) ▲

- ⚠ Warning** Use of controls or adjustments or performance of procedures other than those specified in your printed and on-line documentation may result in hazardous radiation exposure. Do not alter or attempt to remove the laser head from its protective housing or attempt to remove the shields that surround the laser head. Exposure to laser energy or high voltage may result. ▲

- ⚠ Warning** The laser energy from the excitation lasers may be visible or invisible. ▲

- ⚠ Warning** Never stare directly into the laser beam or at its bright reflection. Never tamper with the laser head. Exposure to laser energy or high voltage may result. ▲

- ⚠ Warning** If a cover is loose or distorted or if the interlocks are not operating properly, turn off the system, remove the power cord, and contact your service representative for assistance. ▲
- ⚠ Warning** Avoid skin or eye contact with the radiation produced by the laser beam or any reflections of that beam. If you wish to observe service personnel, wear properly rated laser safety goggles whenever the instrument power is on and the protective covers are open or removed. ▲
- ⚠ Caution** Do not attempt to service the lasers yourself. Please call our office in your area if a laser in your instrument needs service. ▲
- ⚠ Warning** To avoid severe burns from exposure to laser radiation, never put anything except samples in the instrument or against the tip of a fiber optic probe. ▲
- ⚠ Warning** Defeating the safety interlocks on your instrument can expose you to dangerous visible and/or invisible Class IIIb laser radiation at levels up to 500 mW. Do not attempt to defeat the interlocks on the instrument. If a sampling procedure, or some other circumstance, requires the interlocks to be defeated, an interlock defeat kit is available from us. When installed, this defeat kit holds the beam shutters in the open position. ▲

Laser emissions

The accessible radiation levels are below limits defined by the United States Department of Health and Human Services. The laser rating is indicated on the label on the front of the instrument. Using some instrument configurations or accessories (such as the fiber optic accessory), can make your instrument a class IIIb laser product.

Manufacturer's laser information

The following table has information from the laser manufacturer that might be useful if you need to register your instrument.

Laser	Characteristic	Specification
532 nm excitation (green)	type of laser	frequency doubled Nd:YVO4 DPSS
	wavelength	532 nm
	minimum power	15 mW*
	nominal power	24 mW*
	maximum power	30 mW*
	beam diameter	1.9 mm (1/e ²)*
	beam divergence	<8 mrad
	operating voltage	5 VDC
	CDRH classification	Class IIIb
633 nm excitation (red)	type of laser	helium neon (HeNe)
	wavelength	632.8 nm
	minimum power	10 mW*
	nominal power	21 mW*
	maximum power	75 mW*
	beam diameter	1.24 mm (1/e ²)*
	beam divergence	0.68 mrad
	spacing C/2L	165 MHz
	operating voltage	4,400 – 53,000 VDC
CDRH classification	Class IIIb	
780 nm excitation (invisible, low power)	type of laser	externally stabilized diode laser
	wavelength	780 nm
	minimum power	10 mW*
	nominal power	35 mW*
	maximum power	50 mW*
	beam diameter	3.0 mm (1/e ²)*
	beam divergence	1 mrad
	operating voltage	12 VDC
	CDRH classification	Class IIIb
780 nm excitation (invisible, high power)	type of laser	internally stabilized diode laser
	wavelength	780 nm
	minimum power	500 mW*
	nominal power	550 mW*
	maximum power	600 mW*
	operating voltage	12 VDC
CDRH classification	Class IIIb	

* At output of laser head.

Note The laser disclosure notice that accompanied the instrument provides information about the laser(s) included with that instrument. ▲

Laser safety goggles

Exposure to radiation from the excitation laser can result in serious injury and/or blindness. To avoid serious injury, wear laser safety goggles whenever:

- You operate a Class IIIb laser product.
- You enter a laser control area.
- You observe trained personnel servicing the spectrometer.
- You are operating a spectrometer that has the interlocks disabled.
- A fiber optic launcher is installed in the spectrometer.

The following list gives the specifications for protective eyewear for the excitation lasers that are currently available from us:

- If you are using the 780 nm NIR diode laser, wear goggles that have an optical density of 6 and are rated for wavelengths of 780 nm.
- If you are using the 633 nm red laser, the required optical density for laser safety eyewear is 0, so safety goggles are not necessary.
- If you are using the 532 nm green laser, wear goggles that have an optical density of 6 and are rated for wavelengths of 532.

Warning

Always verify that the wavelength(s) listed on your eyewear correspond to the wavelength(s) of the laser being used. ▲

Danger

Before the instrument is used, the laser safety officer should obtain appropriate laser safety eyewear and verify that it is rated for the correct wavelength(s) and optical density for the laser(s) being used. ▲

Disabling the excitation laser

Your instrument has a key switch that allows you to disable the laser. For more information, see the on-line instrument help.

Protective housing

The protective housings on the instrument prevent exposure to laser energy. During normal operation, leave all the covers in place.

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© 2008 Thermo Fisher Scientific Inc. All rights reserved. 本書の全部または一部を弊社の書面による許可なく検索システムに保存したり、または複写、複製したりすることは、 photocopy、写真、磁気その他の記録を含め、さらにこれに限定されないいかなる媒体をもってしても許可されません。